


PURPOSE

To guide fluid resuscitation in the initial management of pediatric trauma.

PROTOCOL

Fluid Resuscitation of the Pediatric Trauma Patient



Notes:

* Rapid Infuser if greater than 25 kg; warm syringe and bolus in less than 5 minutes if less than 25 kg

[△]MTP = Massive Transfusion Protocol

REFERENCES

- Bouglé A, Harrois A, Duranteau J. (2013). Resuscitative strategies in traumatic hemorrhagic shock. *Annals of Intensive Care*, 3(1):1-9.
- Revell, M., Greaves, I. and Porter, K. (2003). Endpoints for Fluid Resuscitation in Hemorrhagic Shock. *Journal of Trauma-Injury Infection & Critical Care*, 54(5):S63-S87.
- LA County & UCLA Trauma Surgery and Critical Care Protocols
- Stainsby, D., MacLennan, S. and Hamilton, P.J. (2000). Management of massive blood loss: A template guideline. *British Journal of Anaesthesia*, 85(3): 487-91
- Schwartzberg, S.D., Bergman, K.S. and Harris, B.H. (1988) A pediatric trauma model of continuous hemorrhage. *Journal of Pediatric Surgery*, 23(7):605–609.